


Patella SuturePlate™ II Fracture Management

Surgical Technique



Introduction

Patella fractures represent approximately 1% of all fractures¹ and present a variety of fracture patterns, which can make osteosynthesis challenging. These fractures are the most common cause of disruption of the extensor mechanism and can result from direct, indirect, or combined trauma.² The amount of comminution, degree of osteoporosis, and the natural forces across the patella create obstacles for reduction and adequate stabilization. Studies have shown significantly higher patellar stability and osteosynthesis rates with locked plates.¹

References

1. Müller EC, Frosch KH. Functional outcomes of revision osteosynthesis after failure of surgical treatment of patellar fractures. *J Knee Surg.* 2021;34(1):80-86. doi:10.1055/s-0039-1692673
2. Pengas IP, Assiotis A, Khan W, Spalding T. Adult native knee extensor mechanism ruptures. *Injury.* 2016;47(10):2065-2070. doi:10.1016/j.injury.2016.06.032

Design and Application

To address various types of fractures, the system offers three different types of plates in different sizes. All plates are made of titanium, and have a low-profile thickness of 1.6 mm with screws sitting flush on the plate. Suture holes in the plate provide soft-tissue reattachment or ligament bracing.

Arrow Plate

Transverse fractures



Star Plate

Comminuted fractures



Star Plate for Pole Fractures

Detached and comminuted distal patella pole




Patient Positioning and Surgical Approach

Patient Positioning

The patient is in supine position with the operated leg prepped and draped.


Surgical Approach

An anterior approach to the knee using a longitudinal midline surgical incision is standard. For simple fracture patterns, the dorsal cortical edges are used to guide reduction of the articular surface. In more complex cases, anterior comminution of the patella may limit cortical reads for reduction. In this scenario, palpation or direct visualization of the articular surface may be required to assure proper reduction. To provide proper access to the articular surface, a lateral arthrotomy is preferred, which preserves the major blood supply of the patella coming inferomedially. For complete visualization of the articular surface, the arthrotomy can be extended to allow a 90-degree eversion of the patella. After fixation of the patella, any traumatic or surgical insult to the retinaculum should be repaired to provide additional strength to the extensor mechanism repair.





Patella SuturePlate™ II Star Plate

Plate Features




The Patella SuturePlate II titanium locking plate is used with 3.0 mm variable-angle locking (VAL) screws. Suture holes in the plate provide soft-tissue reattachment or ligament bracing. In addition to the plate, QuickFix™ 4.0 mm cannulated screws can be used to apply compression through certain fragments.


Radiological Images of the Patella SuturePlate II Star Plate




Star Plate Surgical Technique



Reduce the bone fragments to the anatomic position using large Weber clamps, small tenaculums, or a patella-positioning clamp. A QuickFix™ 4.0 mm cannulated screw can be used to apply compression through certain fragments. Use a drill guide for a 1.35 mm guidewire with trocar tip. Identify the desired screw length using the 4.0 mm cannulated depth gauge. Use the drill guide to overdrill the 1.35 mm guidewire with a cannulated 2.6 mm drill bit.




Use the trial plates with a positioning handle to determine the correct plate size and position. Fluoroscopy of the trial can help define the sizing of the implant. After trialing, open the corresponding sterile implant.




Threaded BB-Taks help to temporarily fix the plate onto the bone. Additionally, locking bending guides can be used to create a more flush fit.

Star Plate Surgical Technique (Cont.)



Screw a drill/depth guide into the locking holes before drilling. Screw length is determined by referencing the laser line on the 2.0 mm drill bit at the scale on the drill/depth guide.

Note: A standard depth gauge can also be used to verify the correct screw length.



Insert and lock the 3.0 mm VAL screws into the patella fracture plate. Fix the locking screws with a T10 hexalobe driver and a driver handle.



Screws must not enter the articulating surface of the patella. Screw heads should sit flush with the plate.

Patella SuturePlate™ II Star Plate for Pole Fractures

Plate Features



The Patella SuturePlate II titanium locking plate is used with 3.0 mm VAL screws. Suture holes in the plate provide soft-tissue reattachment or ligament bracing. In addition to the plate, QuickFix™ 4.0 mm cannulated screws can be used to apply compression through certain fragments. The pole plate can be used in certain fracture patterns in which the distal pole is displaced. Among patellar fractures treated surgically, approximately 20% involve the inferior patellar pole.¹


Radiological Image of the Patella SuturePlate II With Pole Plate




Reference

1. Egol K, Howard D, Monroy A, Crespo A, Tejwani N, Davidovitch R. Patella fracture fixation with suture and wire: you reap what you sew. *Iowa Orthop J.* 2014;34:63-67.

Star Plate for Pole Fractures Surgical Technique




Reduce the bone fragments to the anatomic position using large Weber clamps, small tenaculums, or a patella positioning clamp. A QuickFix™ 4.0 mm cannulated screw can be used to apply compression through certain fragments. Use a drill guide for a 1.35 mm guidewire with trocar tip. Identify the desired screw length using the 4.0 mm cannulated depth gauge. Use the drill guide to overdrill the 1.35 mm guidewire with a cannulated 2.6 mm drill bit.




After identification of the distal pole, perform a double incision of the patellar ligament. Using a positioning handle, position the trial of the patella SuturePlate™ II for pole fractures onto the anterior patella surface, addressing the fracture fragments and including the distal pole. Fluoroscopy of the trial can help define the sizing of the implant. After trialling, open the corresponding sterile implant.

Star Plate for Pole Fractures Surgical Technique (Cont.)




For correct positioning, attach the hook to the distal pole using the double incision of the patellar ligament and place the plate on the anterior cortex of the patella. In addition, the locking bending guide can be used to bend the plate. Threaded BB-Taks may help to temporarily fix the plate onto the bone.

Note: The hooks of the plate can be used to reduce the fracture. The hooks need to be carefully placed at the distal pole. Aggressive impaction of hooks can cause secondary dislocations or additional fractures of the distal pole.



A drill/depth guide is screwed into the locking holes before drilling. Make sure that there is no K-wire or screw in the drilling path. Screw length is determined by referencing the laser line on the 2.0 mm drill bit at the scale on the drill/depth guide.


Note: A standard depth gauge can also be used to verify the correct screw length.



Insert and lock the 3.0 mm VAL screws into the patella pole fracture plate. Fix the locking screws with a T10 hexalobe driver and a driver handle.

Note: Screws must not be placed into the fracture lines.

Star Plate for Pole Fractures Surgical Technique



Screws should not enter the articulating surface of the patella. Screw heads should sit flush with the plate.



Additional Stabilization

In addition to plating the patella, FiberTape® cerclage can also be used for additional stabilization, depending on fracture pattern.

Note: Over-reinforcing the retinaculum by sutures may compromise the blood supply of the patella.

Ordering Information

Implants 3 mm System

Product Description	Item Number
Patella SuturePlate™ II Arrow, sterile	AR-13070A-S
Patella SuturePlate II Star, small, sterile	AR-13070S-S
Patella SuturePlate II Star, medium, sterile	AR-13070M-S
Patella SuturePlate II Star, large, sterile	AR-13070L-S
Patella SuturePlate II Star Plate Pole Fracture, small, sterile	AR-13070S-P-S
Patella SuturePlate II Star Plate Pole Fracture, medium, sterile	AR-13070M-P-S
Patella SuturePlate II Star Plate Pole Fracture, large, sterile	AR-13070L-P-S

Implant Trials

Product Description	Item Number
Patella SuturePlate II Trial Arrow	AR-13070A-T
Patella SuturePlate II Trial Star, small	AR-13070S-T
Patella SuturePlate II Trial Star, medium	AR-13070M-T
Patella SuturePlate II Trial Star, large	AR-13070L-T
Patella SuturePlate II Trial Star, pole fracture, small	AR-13070S-P-T
Patella SuturePlate II Trial Star, pole fracture, medium	AR-13070M-P-T
Patella SuturePlate II Trial Star, pole fracture, large	AR-13070L-P-T

Ordering Information (Cont.)

Screws

Product Description	Item Number
VAL Screw, Ti, 3 mm × 10 mm-40 mm, qty. 4	AR-8933V-10 – 40
Compression Screws	
QuickFix™ Screw, cannulated shaft, cancellous, titanium, 4 mm × 40 mm	AR-8740-40PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 42 mm	AR-8740-42PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 44 mm	AR-8740-44PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 46 mm	AR-8740-46PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 48 mm	AR-8740-48PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 50 mm	AR-8740-50PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 52 mm	AR-8740-52PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 54 mm	AR-8740-54PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 56 mm	AR-8740-56PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 58 mm	AR-8740-58PTS
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 60 mm	AR-8740-60PTS
Optional (not included in the screw caddy)	
QuickFix Screw, cannulated shaft, cancellous, titanium, 4 mm × 28 mm-38 mm	AR-8740-28PTS – 38PTS

Instruments 3 mm System

Product Description	Item Number
Patella Repositioning Clamp	AR-13055
Positioning Handle	AR-14024
Bending Guide, locking, 3 mm	AR-8950-09
TRIM-IT™ Depth Gauge, small	AR-4166
Driver for 3 mm Locking Screws, T10 hexalobe	AR-8944DH
Drill Bit, 2 mm	AR-8944-22
Drill Guide, locking, threaded, 3 mm	AR-8950-07
Handle QC, ratcheting, cannulated	AR-8950RH
Weber Clamp	AR-8943-24
Bone Reduction Forceps, curved, pointed, qty. 2	AR-8943-07

Instruments (for QuickFix Screws)

Product Description	Item Number
Guidewire, with trocar tip, Ø 1.35 mm	AR-8737-01
Depth Guide, cannulated, 4 mm	AR-8737-10
Drill Guide, 2.6 mm/1.35 mm	AR-8943-03
Drill Bit, cannulated, 2.6 mm, qty. 2	AR-8943-02
Driver Shaft, for QuickFix screws, cannulated, T15 hexalobe	AR-8943-12
Holding Sleeve, for 2.7, 3.5 and 4 mm screws	AR-8943-11

Disposables

Product Description	Item Number
K-Wire, 1.6 mm × 150 mm	AR-14016
BB-Tak, threaded	AR-13226T
BB-Tak	AR-13226

Products advertised in this brochure/surgical technique guide may not be available in all countries. For information on availability, please contact Arthrex Customer Service or your local Arthrex representative.



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.

View U.S. patent information at www.arthrex.com/corporate/virtual-patent-marking

© 2021 Arthrex, Inc. All rights reserved. | arthrex.com | LT1-000206-en-US_A