


# Pilon Fusion System

Surgical Technique



---

## Introduction

The Arthrex Pilon Fusion System was designed to treat distal tibia fractures that require not only fracture reduction but also primary ankle arthrodesis. Severe damage to the tibiotalar (TT) joint often results in posttraumatic arthritis, pain, stiffness, and the need for secondary surgeries. The Primary Pilon Fusion System provides another option to address these severe fracture patterns with primary TT arthrodesis of the articular surface to avoid secondary surgery and chronic pain. Anterolateral and posterior approaches, depending on the fracture pattern, allow for fracture management, lengthy bridging techniques, anatomic implants, and fracture-specific locking configurations.

Arthrex offers comprehensive solutions to treat these patients with the ArthroFX™ Large External Fixation System, Ankle Fusion Plating System, Titanium Ankle Fracture Management System, Arthrex Distal Tibia Plating System, and the FibuLock® fibular nail.






## Anterolateral Pilon/Tibiotalar Fusion Plate Design

The anterolateral plate has two rows of distal 3.5 mm locking or 3.5 mm cortical screws to address complex, high-energy pilon fractures. The use of 3.5 mm screws distally allows for a high density of screw fixation and additional options for fracture reduction. The tibial shaft and talus fixation points can use a combination of 4.5 mm locking, 4.5 mm cortical, and 5.5 mm cancellous screws for increased strength across the ankle joint. An oblong slot and K-wire/BB-Tak holes facilitate proper plate placement and provisional fixation.

The anatomic distal contour allows for robust fixation across the ankle joint while providing the maximum number of fixation points for a successful fusion.


- Plate lengths: 5-hole (112 mm), 7-hole (137 mm), 9-hole (163 mm), and 12-hole (201 mm)
- Sterile options: 15-hole (239 mm) and 18-hole (277 mm)




## Anterolateral Primary Pilon Fusion Plate Technique



Use an anterior or anterolateral approach, ensuring the neurovascular structures are protected. The fracture can be addressed with provisional fixation to aid in restoring the articular surface of the distal tibia.




Fracture management and joint preparation are managed accordingly based on the fracture pattern and surgeon preference. The ankle joint can be prepped for a fusion by using the various chisels and curettes offered in the system.




Attach the 4.5 mm percutaneous insertion handle to the most distal shaft screw of the plate. Slide the plate proximally in the previously created submuscular pocket. It is recommended to slide the plate proximally first and then seat it distally into the appropriate position based on anatomic contours.

The percutaneous insertion handle is cannulated to enable provisional fixation with a K-wire if desired.


## Anterolateral Primary Pilon Fusion Plate Technique (Cont.)



Once the plate is positioned, drill with the 3.0 mm drill bit through the 3.0 mm/4.5 mm drill guide. Measure the appropriate screw length with the depth gauge and place a 4.5 mm cortical screw in the oblong slot; conversely, a BB-Tak can be used to provide provisional fixation.




Provisionally fix the talus with a BB-Tak. Drill with the 3.0 mm calibrated drill through the 4.5 mm locking drill guide; the screw length can be measured off the back of the drill guide. Using the T20 driver, insert the appropriate locking or cortical screw.



Use 1.6 mm K-wires to secure the plate to the distal tibia and to secure the fracture fragments. Use a 3.5 mm cortical screw to manipulate the bone toward the plate or a 3.5 mm locking screw to maintain spatial positioning.


Drill with the 2.5 mm calibrated drill through the 3.5 mm locking drill guide. Measure with the drill guide or the depth gauge and implant a 3.5 mm cortical or locking screw with the T15 driver.

## Anterolateral Primary Pilon Fusion Plate Technique (Cont.)




Once the distal fixation is completed add a 4.5 mm locking screw to the kickstand screw hole to provide additional stability across the joint. Using the 4.5 mm locking drill guide drill with the 3.0 mm calibrated drill and measure the screw length off the drill guide.

**Note:** A 4.5 mm nonlocking cortical or 5.5 mm cancellous locking screw may be used if desired.



Distal screw trajectories.



Once the distal fixation is completed add 4.5 mm cortical nonlocking or locking screws proximally as needed.



**Optional:** AlloSync™ demineralized bone graft hydrated with concentrated bone marrow aspirate can be used to augment the ankle arthrodesis procedure.





## Posterior Pilon/Tibiotalar Fusion Plate Design

The posterior plate has two rows of distal 4.5 mm locking or 4.5 mm nonlocking cortical screws to address complex, high-energy pilon fractures while providing stability and strength for a successful fusion. With three points of fixation in the talus, along with the locking kickstand screw across the ankle joint, the posterior plate provides a robust buttress effect for addressing fracture reduction while providing the spatial arrangement and stability needed to maintain talus positioning for tibiotalar arthrodesis.

An oblong slot and K-wire/BB-Tak holes facilitate proper plate placement and provisional fixation. The anatomic distal contour allows for robust fixation across the ankle joint while providing the maximal number of fixation points.


- Plate lengths: 5-hole (119 mm), 7-hole (144 mm), 9-hole (170 mm), 12-hole (208 mm)
- Sterile options: 15-hole (246 mm)



## Posterior Pilon/Tibiotalar Fusion Plate Technique




Use a posterolateral or posteromedial approach, dictated by fracture pattern and surgeon preference, while carefully protecting tendons and neurovascular structures. Obtain provisional fracture reduction and prepare the ankle joint surfaces for fusion.




Attach the 4.5 mm percutaneous insertion handle to the most distal shaft screw of the plate. Slide the plate proximally past the ideal insertion point and then move it distally until it seats in the appropriate position on the distal tibia poster facet of the talus.

The percutaneous insertion handle is cannulated to enable provisional fixation with a K-wire if desired.


## Posterior Pilon/Tibiotalar Fusion Plate Technique (Cont.)




Place a BB-Tak in the oblong hole for provisional fixation and begin fixation with 4.5 mm cortical screws. Drill with the 3.0 mm drill bit through the 3.0/4.5 mm drill guide.



Measure the length with the depth gauge and insert the desired 4.5 mm screw. Remove the BB-Tak.




Add additional 4.5 mm screws as needed. Strategically placed screws may help with manipulations of the fracture fragment position.




Secure the talus by placing a BB-Tak in the medial or lateral talar screw hole. A 4.5 mm cortical screw can help decrease any gaps between the bone and the plate, but may affect alignment. Be sure to maintain the desired foot and fracture alignment while drilling. The use of 4.5 mm locking screws will help maintain fracture alignment. Drill with the 3.0 mm drill bit and the 3.0 mm/4.5 mm drill guide.

## Posterior Pilon/Tibiotalar Fusion Plate Technique (Cont.)




Once the distal fixation is completed, add a 4.5 mm locking screw to the kickstand screw hole to provide additional stability across the joint. Use the 4.5 mm locking drill guide drill with the 3.0 mm calibrated drill and measure the screw length off the drill guide.




Distal screw trajectories.

**Note:** The locking screw trajectories of the kickstand screw and the central talar screw will intersect. The surgeon must be mindful during drilling and select screws of the appropriate length to avoid interference.






Once the distal fixation is completed, add a 4.5 mm cortical or locking screw as needed.



**Optional:** AlloSync™ demineralized bone graft hydrated with concentrated bone marrow aspirate can be used to augment the ankle arthrodesis procedure.



## Tray Overview



# Supporting Products



## Arthrex Ankle Fusion Plating System

The titanium Ankle Fusion Plating System provides a complete solution for ankle fusion management with a comprehensive offering of anatomy-specific plates available for either tibiotalar or tibiotalocalcaneal arthrodesis. A variety of screw options, including locking, nonlocking, cortical, cancellous, and hybrid designs, are provided to address all fixation needs. Specific instrumentation designed to help gain access to and prepare the fusion sites is included in the set for completeness. The Ankle Fusion Plating System was designed to provide the solution to your ankle fusion fixation needs.


### System Features

- Anatomically designed for use with three surgical approaches: anterior, lateral, posterior
- Four compression modes available in system
  - Anatomic compression hole
  - Oblong compression hole
  - Mini joint compressor/distractor
  - 6.7 mm cannulated lag screws or 7.0 mm XL Compression FT screws


- Comprehensive instrumentation for joint preparation, distraction, and compression and assistance with optimal fixation
- Maximum fixation points within each plate




Anterior Tibiotalar Plate  
AR-8970AR



Lateral Tibiotalar Plate  
AR-8970TT



Lateral Tibiototalcaneal  
Plate – AR-8970TTC



Posterior Tibiototalcaneal  
Plate – AR-8970PR




### Cannulated Screws

The versatility of the Ankle Fusion System provides a complete solution for treating ankle arthritis in one comprehensive instrument case. The instrument set can be configured to house either 6.7 mm cannulated lag screws or 7.0 mm XL Compression FT screws for percutaneous compression across the arthrodesis site.




6.7 mm LPS Screw (18 mm thread)



7.0 mm XL Compression FT Screw

- **Low-profile head** – 1 mm shorter than a traditional 6.5 mm AO screw, while still using a 3.5 mm hex driver.
- **Increased pull-out** – 30% better than a standard 6.5 mm AO screw.<sup>1</sup>
- **Deeper threads** – Using a 2.4 mm guide pin allows the threads to be deeper than a standard AO screw.
- **Self-drilling/tapping** – Speeds up the insertion process.
- **Headless design** – Minimal risk of impingement or soft-tissue irritation.
- **Fully threaded compression** – Variable-stepped thread pitch and tapered proximal profile work together to compress bone fragments with the purchase of a fully threaded screw.
- **Self-drilling/tapping** – Helical relief flutes assist in bone removal to reduce insertion torque.



Osteotome Angled Up, 5.5 mm

Osteotome Angled Up, 12 mm

Osteotome Straight, 5.5 mm

Osteotome Straight, 12 mm

Cup Curette, straight, 6 mm

Cup Curette, curved, 6 mm

### Joint Preparation

Straight and curved curettes and osteotomes have been added to the ankle fusion tray to help with the removal of cartilage from the ankle and subtalar joints. These instruments are appropriately designed for the ankle and come standard in each ankle fusion tray, simplifying joint preparation in the OR setting.

### Mini Joint Compressor/Distractor


Adaptable for distraction and compression of arthrodesis sites, this unique device facilitates joint preparation and allows for excellent compression prior to definitive fixation. The device uses 1.6 mm or 2.4 mm guidewires, or 3 mm traction screws, which are included in the system.




### Reference

1. Robert KQ, Chandler R, Baratta RV, et. al. The effect of divergent screw placement on the initial strength of plate-to-one fixation. *J Trauma*. 2003;55(6):1139-1144. doi:10.1097/01.TA.0000031103.15337.CA







**Level 1**  
Instrument Case  
Implant Tray and  
Instruments  
AR-8970C-01



**Level 2**  
Interchangeable  
Cannulated Screw  
Instrument Tray




**Level 2**  
**7.0 mm Screw Instrumentation**  
■ AR-8970C-70 – 7.0 mm Tray  
■ AR-8970C-SC-70 – 7.0 mm  
Screw Caddy Insert



**Level 2**  
**6.7 mm Screw Instrumentation**  
■ AR-8970C-67 – 6.7 mm Tray  
■ AR-8970C-SC-67 – 6.7 mm  
Screw Caddy Insert








**Level 3**  
Auxiliary  
Instruments and  
Screw Caddy



**Screw Caddy**  
Interchangeable  
Cannulated Screws



**Screw Caddy Inserts**  
■ AR-8970C-SC-70  
■ AR-8970C-SC-67

# Supporting Products

## Biologic Options

### Angel® Concentrated Platelet-Rich Plasma (cPRP) System


Technology is what sets the Angel system apart from the competition. The Angel system is the only one to provide PRP concentrate from bone marrow aspirate (BMA) with adjustable cellular levels. Bone marrow is a rich source of platelets and nucleated and progenitor cells. Customization of cellular levels is necessary to reduce the number of neutrophils in bone marrow concentrate (BMC), which can be detrimental to bone healing.

#### Features and Benefits:


- Proprietary platelet sensor system
- Adjustable platelet concentrations
- Adjustable WBC concentrations
- Flexible processing volume 40 mL to 180 mL
- Each processing kit can process 3 cycles up to 180 mL, on the same patient
- Programmable – can store up to 30 custom processing protocols
- Closed system: delivers PRP, platelet-poor plasma, and RBCs into separate, sterile compartments



In vitro culture expansion of progenitor cells<sup>1</sup>



48 hours




96 hours

Angel cPRP System	Platelet Concentration (K/ $\mu$ L)	Nucleated Cell Concentration (K/ $\mu$ L)	Hematopoietic Cell Concentration (K/ $\mu$ L)	Total Neutrophils ( $\times 10^6$ )
BMA	87.7 $\pm$ 6.4	24.5 $\pm$ 15.6	0.002 $\pm$ 0.001	612.1
BMC	787.0 $\pm$ 317.6	240.5 $\pm$ 186.6	0.081 $\pm$ 0.056	132.9
Increase Above Baseline	$\sim 9\times$	$\sim 10\times$	$\sim 33\times$	80% $\downarrow$

Data from Arthrex, Inc. Data on file (APT-02569). Naples, FL; 2018.

### AlloSync™ Pure Demineralized Bone Matrix

AlloSync Pure demineralized bone matrix is derived from 100% human allograft bone with no extrinsic carriers. AlloSync Pure bone matrix resists irrigation and can be used in a fluid environment. The clinician can control the handling properties of AlloSync Pure bone matrix, which includes decreasing the viscosity for injectable applications or increasing the viscosity for open procedures. The proprietary rice-shape fiber technology used to process AlloSync Pure bone matrix increases the osteoinduction and osteoconductive surface area to accelerate cellular ingrowth.<sup>2</sup>




#### References

1. Arthrex, Inc. Data on file (APT-05220). Naples, FL; 2021.
2. Martin GJ, Boden SD, Titus L, Scarborough NL. New formulations of demineralized bone matrix as a more effective graft alternative in experimental posterolateral lumbar spine arthrodesis. *Spine*. 1999;24(7):637-645. doi:10.1097/00007632-199904010-00000

## Supporting Products – Biologic Options

### AlloSync™ Demineralized Bone

Demineralized cancellous sponges and cortical fibers are optimal for combination with blood products such as concentrated BMA. When combined with bone marrow concentrate (BMC), AlloSync demineralized bone grafts provide the necessary components for bone formation: cell, signal, and scaffold.<sup>1</sup>



## Ordering Information

### Arthrex Angel® System

Product Description	Item Number
Arthrex Angel System	ABS-10060
Angel Bone Marrow Processing Kit	ABS-10062
Angel Blood Access Kit	ABS-10067

### AlloSync Cancellous Sponges

Product Description	Item Number
Cube, 8 mm × 8 mm × 8 mm	ABS-2005-01
Cube, 10 mm × 10 mm × 10 mm	ABS-2005-02
Cube, 12 mm × 12 mm × 12 mm	ABS-2005-03
Strip, 10 mm × 10 mm × 3 mm	ABS-2006-01
Strip, 15 mm × 40 mm × 3 mm	ABS-2006-02
Strip, 26 mm × 19 mm × 7 mm	ABS-2006-03
Strip, 10 mm × 20 mm × 7 mm	ABS-2006-04
Chips (1 mm-4 mm), 1.0 cc	ABS-2007-01
Chips (1 mm-4 mm), 2.5 cc	ABS-2007-02
Chips (1 mm-4 mm), 5 cc	ABS-2007-03

### AlloSync Cortical Fibers

Product Description	Item Number
Fibers, 1.0 cc	ABS-2008-01
Fibers, 2.5 cc	ABS-2008-02
Fibers, 5 cc	ABS-2008-03
Fibers, 10 cc	ABS-2008-04

### AlloSync Pure

Product Description	Item Number
Pure, 1.0 cc	ABS-2010-01
Pure, 2.5 cc	ABS-2010-02
Pure, 5 cc	ABS-2010-05
Pure, 10 cc	ABS-2010-10

### Reference

1. Kay JF, Khaliq S, Neubauer P. Effective design of bone graft materials using osteoinductive and osteoconductive components. American Association of Tissue Banks. <https://www.aatb.org/sites/default/files/2003Abstract13.pdf>. Accessed January 23, 2018






## Ordering Information

### Primary Pilon Fusion System (AR-8963FS)

Product Description	Item Number
Screwdriver, T15 hexalobe	AR-8943-10
Drill Guide, 3.5/2.5 mm	AR-8943-14
Drill Guide, locking, 3.5 mm	AR-8943-43
Countersink, 3.5/4.0 mm	AR-8950-03
Depth Device, 3.5/4.0 mm	AR-8963-13
Chisel, straight, 14 mm	AR-8963-14
Chisel, curved, 14 mm	AR-8963-15
Chisel, curved, 7mm	AR-8963-16
Chisel, straight, 7 mm	AR-8963-17
Cup Curette, 9 mm	AR-8963-18
Lag Drill Guide, 3.5 mm/2.5 mm	AR-8963-21
Percutaneous Insertion Handle, 4.5 mm	AR-8963-24
Driver, T15 hexalobe, 6 in, AO, qty. 2	AR-8963-25
Driver, T15 hexalobe, 6 in, straight, AO, qty. 2	AR-8963-26
Primary Pilon Fusion Case	AR-8963C-03
<b>Plates</b>	
Anterolateral Pilon Fusion Plate, 5h, left	AR-8963ALF-05
Anterolateral Pilon Fusion Plate, 7h, left	AR-8963ALF-07
Anterolateral Pilon Fusion Plate, 9h, left	AR-8963ALF-09
Anterolateral Pilon Fusion Plate, 12h, left	AR-8963ALF-12
Anterolateral Pilon Fusion Plate, 5h, right	AR-8963ARF-05
Anterolateral Pilon Fusion Plate, 7h, right	AR-8963ARF-07
Anterolateral Pilon Fusion Plate, 9h, right	AR-8963ARF-09
Anterolateral Pilon Fusion Plate, 12h, right	AR-8963ARF-12
Posterior Pilon Fusion Plate, 5h	AR-8963PF-05
Posterior Pilon Fusion Plate, 7h	AR-8963PF-07
Posterior Pilon Fusion Plate, 9h	AR-8963PF-09
Posterior Pilon Fusion Plate, 12h	AR-8963PF-12

Product Description	Item Number
<b>Sterile Plates</b>	
Anterolateral Pilon Fusion Plate, 15h, left	AR-8963ALF-15S
Anterolateral Pilon Fusion Plate, 18h, left	AR-8963ALF-18S
Anterolateral Pilon Fusion Plate, 15h, right	AR-8963ARF-15S
Anterolateral Pilon Fusion Plate, 18h, right	AR-8963ARF-18S
Posterior Pilon Fusion Plate, 15h	AR-8963PF-15S
<b>3.5 mm Screws, Low-Profile, Ti</b>	
3.5 mm × 20 mm-60 mm (2 mm increments)	AR-8935-20-60
3.5 mm × 65 mm-80 mm (5 mm increments)	AR-8935-65-80
3.5 mm × 90 mm-120 mm (10 mm increments)	AR-8935-90-120
<b>3.5 mm Screws, Low-Profile, Ti, Locking</b>	
3.5 mm × 20 mm-50 mm (2 mm increments)	AR-8935L-20-50
3.5 mm × 55 mm-60 mm (2 mm increments)	AR-8935L-55-60
<b>Disposables</b>	
Drill Bit, 3.5 mm, qty. 2	AR-4160-35
Drill Bit, 2.5 mm, qty. 2	AR-8963-19
Guidewire, drill tip, 3.0 mm, qty. 4	AR-8963-20
Drill Bit, calibrated, long, 3.0 mm	AR-8970-30L
Bone Tap	AR-8963-23





## Ordering Information (Cont.)

### Ankle Fusion Plating System, 7.0 mm Set (AR-8970S-70)

Product Description	Item Number
Perc Drill Guide, Compression FT	AR-8750-02
Drill Guide, threaded, locking, 4.5 mm, qty. 2	AR-8970-01
Drill Guide, 3 mm/4.5 mm	AR-8970-02
Drill Guide, 3 mm/5.5 mm	AR-8970-05
Depth Measuring Device, long, 4.5 mm/5.5 mm	AR-8970-07L
Depth Device, cannulated screws	AR-8750-01
Drive Shaft, T20 hexalobe, qty. 2	AR-8970-03
Driver, T20 hexalobe, straight	AR-8970-04
Driver, T20 hexalobe, straight, AO, qty.2	AR-8970-08
Driver, T25 hexalobe, ISO, cannulated, qty. 2	AR-8770-01
Driver, T25 hexalobe, ISO, solid	AR-8770-04
Ratcheting Handle, cannulated, large AO handle, QC	AR-8970RH
Mini Joint Distractor/Compressor	AR-8970JD
Axial Handle, trilobe QC, ratcheting	AR-8770RH
Soft Tissue Protector, 2.4 mm	AR-8770-06
Bone Reduction Forceps, qty. 2	AR-8943-07
Hohmann Retractor, 9.5 in, 17 mm pointed, qty. 2	AR-9260-34
Cup Curette, straight shaft, 6 mm	AR-8970-11
Cup Curette, curved shaft, 6 mm	AR-8970-12
Cobb Elevator, 9 mm	AR-8640
Small Joint Osteotome Angled Up, 0.217 in (5.5 mm) w/ handle	AR-8650-08
Small Joint Osteotome Straight, 0.217 in (5.5 mm) w/ handle	AR-8650-09
Small Joint Osteotome Angled Up, 0.472 in (12 mm) w/ handle	AR-8970-13
Small Joint Osteotome Straight, 0.472 in (12 mm) w/ handle	AR-8970-14
Screw Holding Forceps	AR-8941F
Guidewire Sleeve Insert, 1.6 mm	AR-8970-06
Ankle Fusion Plating System Instrument Case	AR-8970C-01
Ankle Fusion Instrument Case, 7.0 mm tray	AR-8970C-70
Ankle Fusion Caddy, 7.0 mm insert	AR-8970C-SC-70

### Disposables for AR-8970S-70 (not included in set, order separately)

Product Description	Item Number
Drill Bit, calibrated, long, 3 mm	AR-8970-30L
Drill Bit, cannulated, long, 3 mm	AR-8970-30CL
Drill Bit, cannulated, long, 4.5 mm	AR-8970-45CL
Drill Bit, long, 4.5 mm	AR-8970-45L
Drill Bit, long, 5.5 mm	AR-8970-55L
Drill Bit, cannulated, long, 5.5 mm	AR-8970-55CL
Drill Bit, cannulated, 5.0 mm	AR-8770-02
Profile Drill, X-large, 7.0 mm	AR-8770-03
BB-Tak, large	AR-8970-09
BB-Tak, large, threaded	AR-8970-09T
Traction Screw, 20 mm	AR-8950JD-2

### Disposables for AR-8970S-70 (not included in set, order separately)

Product Description	Item Number
Traction Post, threaded, 4.5 mm	AR-8970JD-45S
Guidewire w/ Trocar Tip, .095 in (2.4 mm) × 9.25 in	AR-8770K
Guidewire w/ Trocar Tip, threaded, .094 in (2.4 mm) × 9.25 in	AR-8770KT
Guidewire w/ Trocar Tip, .062 in (1.6 mm) × 7 in	AR-8941-7

### Plates for 6.7 mm/7.0 mm Sets (order separately)

Product Description	Item Number
Anterior Plate, 3H, left	AR-8970AL
Anterior Plate, 4H, left	AR-8970AL-04
Anterior Plate, 5H, left	AR-8970AL-05
Anterior Plate, 6H, left	AR-8970AL-06
Anterior Plate, 3H, right	AR-8970AR
Anterior Plate, 4H, right	AR-8970AR-04
Anterior Plate, 5H, right	AR-8970AR-05
Anterior Plate, 6H, right	AR-8970AR-06
Anterior Plate, short	AR-8970AS-03
Anterior Plate, minimally invasive	AR-8970MA
Lateral Tibiotalar Plate, 3H	AR-8970TT
Lateral Tibiotalar Plate, 4H	AR-8970TT-04
Lateral Tibiotalar Plate, 5H	AR-8970TT-05
Lateral Tibiotalar Plate, 6H	AR-8970TT-06
Lateral Tibiotalocalcaneal Plate, 3H	AR-8970TTC
Lateral Tibiotalocalcaneal Plate, 4H	AR-8970TTC-04
Lateral Tibiotalocalcaneal Plate, 5H	AR-8970TTC-05
Lateral Tibiotalocalcaneal Plate, 6H	AR-8970TTC-06
Posterior Tibiotalocalcaneal Plate, left	AR-8970PL
Posterior Tibiotalocalcaneal Plate, right	AR-8970PR

### Low Profile Screws, 4.5 mm/5.5 mm Screws

Product Description	Item Number
Low Profile Locking Screws	
4.5 mm × 18 mm-50 mm (2 mm increments)	AR-8545L-18-50
4.5 mm × 55 mm-75 mm (5 mm increments)	AR-8545L-55-75
Low Profile Screws	
4.5 mm × 18 mm-50 mm (2 mm increments)	AR-8545-18-50
4.5 mm × 55 mm-100 mm (5 mm increments)	AR-8545-55-100
Low Profile Screws, cancellous	
5.5 mm × 20 mm-100 mm (5 mm increments)	AR-8555-20-100

### 7.0 XL Compression FT screws

Product Description	Item Number
7.0 XL Compression FT Screws, cannulated, Ti, fully threaded, 40 mm-100 mm (5 mm increments)	AR-8770-40H-100H

## Ordering Information (Cont.)

### Ankle Fusion Plating System, 6.7 mm Set (AR-8970S-67)


Product Description	Item Number
Drill Guide, threaded, locking, 4.5 mm, qty. 2	AR-8970-01
Drill Guide, 3 mm/4.5 mm	AR-8970-02
Drill Guide, 4 mm/6.7 mm	AR-8967G
Drill Guide, 3 mm/5.5 mm	AR-8970-05
Depth Measuring Device, long, 4.5 mm/5.5 mm	AR-8970-07L
Depth Device, cannulated, for 6.7 mm screws	AR-8967DG
Depth Device, large	AR-4167
Drive Shaft, T20 hexalobe, qty. 2	AR-8970-03
Driver, cannulated, 3.5 mm hex, qty. 2	AR-8967D
Driver, T20 hexalobe, straight	AR-8970-04
Driver, T20 hexalobe, straight, AO, qty.2	AR-8970-08
Ratcheting Handle, cannulated, large AO handle, QC	AR-8970RH
Mini Joint Distractor/Compressor	AR-8970JD
Screwdriver Handle, ratcheting	AR-1999
Bone Reduction Forceps, qty. 2	AR-8943-07
Hohmann Retractor, 9.5 in, 17 mm pointed, qty. 2	AR-9260-34
Hudson Adapter	AR-1416
Cup Curette, straight shaft, 6 mm	AR-8970-11
Cup Curette, curved shaft, 6 mm	AR-8970-12
Cobb Elevator, 9 mm	AR-8640
Screw Holding Forceps	AR-8941F
Countersink, fixed handle, cannulated, 6.7 mm	AR-8967CSF
Guidewire Sleeve Insert, 1.6 mm	AR-8970-06
Small Joint Osteotome Angled Up, 0.217 in (5.5 mm) with handle	AR-8650-08
Small Joint Osteotome Straight, 0.217 in (5.5 mm) with handle	AR-8650-09
Small Joint Osteotome Angled Up, 0.472 in (12 mm) with handle	AR-8970-13
Small Joint Osteotome Straight, 0.472 in (12 mm) with handle	AR-8970-14
Ankle Fusion Instrument Case	AR-8970C-01
Ankle Fusion Instrument Case, 6.7 mm tray	AR-8970C-67
Ankle Fusion Caddy, 6.7 mm insert	AR-8970C-SC-67

### Disposables for AR-8970S-67 (order separately)

Product Description	Item Number
BB-Tak, large	AR-8970-09
BB-Tak, large, threaded	AR-8970-09T
Guidewire w/ Trocar Tip, nonthreaded, 0.094 in (2.4 mm) × 8 in, qty. 6	AR-8967K
Guidewire w/ Trocar Tip, threaded, 0.094 in (2.4 mm) × 8 in, qty. 6	AR-8967KT
Guidewire w/ Trocar Tip, nonthreaded, 0.094 in (2.4 mm) × 12 in, qty. 6	AR-8967K-12
Guidewire w/ Trocar Tip, threaded, 0.094 in (2.4 mm) × 12 in, qty. 6	AR-8967KT-12
Guidewire w/ Trocar Tip, .062 in (1.6 mm) × 7 in, qty. 6	AR-8941-7
Washer, Ti, 13 mm	AR-8967W
Traction Screw, 20 mm	AR-8950JD-2
Drill Bit, cannulated, 4.0 mm	AR-8970-40C
Drill Bit, calibrated, long, 3 mm	AR-8970-30L
Drill Bit, cannulated, long, 3 mm	AR-8970-30CL
Drill Bit, cannulated, long, 4.5 mm	AR-8970-45CL
Drill Bit, long, 4.5 mm	AR-8970-45L
Drill Bit, long, 5.5 mm	AR-8970-55L
Drill Bit, cannulated, long, 5.5 mm	AR-8970-55CL

### Cannulated Lag Screws

Product Description	Item Number
Low Profile Screws, cannulated, partially threaded, 6.7 mm × 40 mm-100 mm, 18 mm length (5 mm increments)	AR-8967-1840-18100



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience, and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level and/or outcomes.

View U.S. patent information at [www.arthrex.com/corporate/virtual-patent-marking](http://www.arthrex.com/corporate/virtual-patent-marking)

**arthrex.com**