




Double Row Rotator Cuff Repair
using the Bio-Corkscrew® FT

Surgical Technique



Double Row Rotator Cuff Repair

Introduction

In the progression of arthroscopic treatment of rotator cuff tears, the double row arthroscopic rotator cuff repair was developed. This advanced procedure is intended to help reestablish the normal footprint of the rotator cuff, enhance mechanical integrity, and improve healing with better clinical outcomes over single row rotator cuff repair.

This repair can be completed using the 4.5, 5.5, or 6.5 mm Bio-Corkscrew FT Suture Anchors. The fully threaded design of these suture anchors increases pull-out strength and reduces suture "pull-back" in soft bone by engaging both cortical and cancellous bone.

Patient Positioning

The patient may be positioned in the beach chair position using the Beach Chair Lateral Traction Device or in a lateral decubitus position using the 3-Point Shoulder Distraction System. Access to the subacromial space is facilitated with a variety of clear cannulas.



Rotator Cuff Tear Assessment: Margin Convergence Repair




Using a KingFisher[®] Suture Retriever/Tissue Grasper or Rotator Cuff Grasper the mobility of the tear is assessed to determine whether a U or L-shaped component exists. In the case of large tears extending to the superior aspect of the glenoid, irrespective of shape, margin convergence suturing is performed in the following manner to reduce volume and strain on the repair:

Via anterior/posterior portal or percutaneous portals consider:

1. BirdBeak[®] to BirdBeak suture hand-off
2. SutureLasso[™]/FiberStick[™] hand-off to BirdBeak/Penetrator[™]
3. Micro SutureLasso/FiberStick to BirdBeak/Penetrator hand-off
4. Scorpion[™] Suture Passer

Soft tissue releases may be necessary in massive rotator cuff tears. These releases can be performed using Tissue Elevators or straight/curved arthroscopic scissors. Refer to DVD-1069 - Complex Arthroscopic Rotator Cuff Repairs (web only)

Double Row Rotator Cuff Repair




Medial Row Anchor Placement


After assessing the width of the rotator cuff footprint, the most medial row suture anchors are placed adjacent to the articular margin of the humerus.

Patient age and bone quality determine the selection of 4.5, 5.5 or 6.5 mm anchors to secure the medial row of the rotator cuff repair.

Pilot hole preparation in the 45° “deadman” angle with the Bio-Corkscrew FT Punch and optional Tap will assist in determining the most suitable anchor. Generally the harder the bone the smaller size anchor that may be used.



Anchors are placed to assure full contact of the detached tendon along the medial footprint of the greater tuberosity.



Medial Row Suture Passing

With smaller tears, the sutures of the medial anchor may be passed and tied in horizontal mattress configuration through the anterior and posterior portion of the tear.


The modified Neviaser portal and other medially based percutaneous portals are considered to pass horizontal mattress sutures just lateral to the musculotendinous junction of the tendon using a Banana SutureLasso, Micro SutureLasso, Banana BirdBeak Evolution or Penetrator for the medial row.

The Scorpion or NeedlePunch II Suture Passers may be used from the lateral portal and will provide the following depth of suture passage through the tissue:

Scorpion, 16 or 20 mm
NeedlePunch II, 10 or 16 mm




Note: Sutures passed through the medial row are preferably stored and tied after the lateral row is completed and the cuff is tensioned to the lateral margin of the footprint.



Lateral Row Anchor Placement

Using the 45° “deadman” angle for optimal anchor insertion, pilot hole preparation for the Bio-Corkscrew FT 4.5, 5.5 mm or 6.5 mm is carried out at the far lateral portion of the rotator cuff footprint and greater tuberosity.



anchors are placed in a linear fashion from anterior to posterior to assist in organizing sutures. Opposite color FiberWire® sutures are alternately passed through tissue and tied in sequence from posterior to anterior.



Lateral Row Suture Passing

Using the lateral portal, vertical mattress sutures are passed 5 mm apart and 10 mm from the tendon edge using primarily the Scorpion Suture Passer or NeedlePunch II.


Alternatively, the anterior and posterior portals are used to pass sutures in the lateral tendon edge with a curved SutureLasso.

(See the Scorpion technique guide video, DVD-1074)

The lateral edge of the tendon is secured to the bone using low profile sliding knots. The medial row sutures are then secured in a similar fashion.

Knot Tying


For double row repairs, it is preferred that reduction and securing the repair to the lateral edge of the rotator cuff attachment site is carried out prior to securing the medial based sutures passed in horizontal mattress fashion. The Weston Knot is a low profile sliding knot and is performed as shown below. In the event of sutures not sliding easily in more complex suture passing situations, alternating half-hitches with the 6th Finger or Single-Hole Knot Pusher will yield secure and low profile knots with #2 FiberWire or TigerWire® sutures.



Crescent-Shape Tears




1 The tear pattern and mobility are assessed (as previously mentioned). Place the required number of medial row anchors, pass the sutures in a mattress stitch configuration and tie the medial row.




2 Pass the required number of lateral row anchors, pass the sutures in a simple stitch configuration and tie the lateral row to complete the repair.

L or Reverse L-Shape Tears



1 The tear pattern and mobility of the tear are assessed (as previously mentioned). The first side-to-side stitch is placed at the lateral edge of the tear, thus may be incorporated within a Bio-Corkscrew FT anchor. Pass the remaining required stitches to complete the margin convergence.



2 The medial and lateral rows of anchors are placed and sutures passed as previously discussed to complete the repair.

Product Information

Bio-Corkscrew FT, 4.5 x 15 mm, w/two #2 FiberWire
 Bio-Corkscrew FT, 4.5 x 15 mm, w/two #2 TigerTail
 Bio-Corkscrew FT, 4.5 x 15 mm, w/Needles and two #2 FiberWire
 Punch for 4.5 mm Corkscrew FT
 Punch/Tap for 4.5 mm Corkscrew FT
 Disposable Punch for 4.5 mm Corkscrew FT

Bio-Corkscrew FT, 5.5 mm x 15 mm, w/two #2 FiberWire
 Bio-Corkscrew FT, 5.5 mm x 15 mm, w/two #2 FiberWire & Needles
 Bio-Corkscrew FT, 5.5 mm x 15 mm, w/two #2 TigerTail
 Bio-Corkscrew FT, 5.5 mm x 15 mm, w/three #2 FiberWire
 BioComposite Corkscrew FT, 5.5 mm x 15 mm, w/two #2 FiberWire
 BioComposite Corkscrew FT, 5.5 mm x 15 mm, w/two #2 TigerTail
 BioComposite Corkscrew FT, 4.5 mm x 15 mm, w/two #2 FiberWire
 Bio-Corkscrew FT w/four NeedlePunch Needles, 5.5 mm x 15 mm, w/two #2 FiberWire
 Bio-Corkscrew FT, 5.5 mm x 15 mm, w/three #2 FiberWire
 Bio-Corkscrew FT Punch for 5.5 mm
 Bio-Corkscrew FT Punch, disposable for 5.5 mm
 Punch/Tap for 5.5 mm Bio-Corkscrew FT

Bio-Corkscrew FT, 6.5 mm x 15 mm, w/two #2 FiberWire

Instrumentation for Suture Passing:

Scorpion Suture Passer, 16 mm
 Scorpion Suture Passer, 20 mm
 Humpback Scorpion Suture Passer, 16 mm
 Scorpion Needle
 NeedlePunch II, 10 mm
 NeedlePunch II, 16 mm
 Suture Shuttle (for NeedlePunch II)

BirdBeak Evolution, 15° up curve
 Banana BirdBeak Evolution, 22° up curve
 Penetrator Suture Retriever, straight

Banana SutureLasso
 Micro SutureLasso

KingFisher Suture Retriever/Tissue Grasper
 Rotator Cuff Grasper
 Crystal Cannula, 5.75 mm I.D. x 7 cm

Media:

Complex Arthroscopic Rotator Cuff Repairs (web only)
 Arthroscopic Rotator Cuff Repair featuring the Scorpion Suture Passer

AR-1927BF-45
 AR-1927BFT-45
 AR-1927BNF-45
 AR-1922P
 AR-1927PTB-45
 AR-1922PBS

AR-1927BF
 AR-1927BNF
 AR-1927BFT
 AR-1927BF-3
 AR-1927BCF
 AR-1927BCFT
 AR-1927BFT-45
 AR-1927BNP4
 AR-1927BF-3
 AR-1927PB
 AR-1927PBS
 AR-1927CTB

AR-1927BF-65

AR-13990
 AR-13992
 AR-13993
 AR-13990N
 AR-13981S
 AR-13982
 AR-7224

AR-11881E
 AR-11892E
 AR-2167ST

AR-4065B
 AR-8702

AR-13970SR
 AR-13960
 AR-6560

DVD-1069
 DVD-1074



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's Directions For Use.

U.S. PATENT NOS. 5,964,783; 6,074,403; 6,517,552; 6,716,234; 6,994,719; 7,029,490 and PATENT PENDING

© 2009, Arthrex Inc. All rights reserved. LT0215C