




Mini Open Latarjet

Surgical Technique




Mini Open Latarjet Technique




A 5 cm skin incision is made starting at the tip of the coracoid process and extending inferiorly, through the deltopectoral approach. Mayo scissors are used to clear the superior aspect of the coracoid process and a Hohmann retractor is placed over the top of the coracoid process.


- Incise clavicopectoral fascia
- Dissect CA ligament
- Release the pectoralis minor from the coracoid



- Dissect fibrofatty tissue from inferior coracoid to its base
- Dissect up to CC ligaments
- Use 90° angled saw blade or osteotome to resect the coracoid




- Prepare the side of the coracoid that will face the glenoid
- Drill two parallel holes in the coracoid using the coracoid drill guide and a 2.75mm drill




Wedged Profile Plate

- The wedged profile improves the match between the coracoid bone block and the glenoid
- Four spikes for more stability
- Two screw holes distribute the load evenly and avoid fracture of the bone block


NOTE: The thicker part of the plate has to be positioned medial to the glenoid neck.




- Choose the cannulated temporary compression device and fix one side of the plate to the coracoid
- Define the superior border of the subscapularis
- Split the subscapularis and open the capsule to expose the anterior glenoid
- Use the Gelpi retractor to keep the joint open
- Use the Fukuda, Swan and Blade retractors to ensure the access to the anteroinferior part of the glenoid




- Prepare the glenoid neck with a burr




- Position the coracoid on the prepared glenoid with the temporary compression device and drill one guidewire through the free hole of the plate into the glenoid
- Drill a second guidewire through the cannulation of the temporary compression device





- Measure the appropriate length of the screw needed
- Leave the temporary compression device in position
- Insert the first screw over the guidewire into the plate/bone construct



- Make sure the bone graft is positioned properly on the glenoid rim
- Remove the temporary compression device and apply the second screw over the guidewire
- The plate ensures a proper contact area between coracoid and glenoid



- Conjoined tendons provide extra joint stability
- Labrum can be refixed to the glenoid using suture anchors which should be placed next to the graft




- The CA ligament stump is sutured to the capsule
- Layer-by-layer closure of the wound


Alternative Technique



- Dissect fibrofatty tissue from inferior coracoid to its base
- Dissect up to CC ligaments
- Use osteotome or 90° angled saw blade to resect the coracoid




- Prepare the side of the coracoid that will face the glenoid
- Drill two parallel holes in the coracoid using the coracoid drill guide and a 4 mm drill




- Select the appropriate glenoid drill guide (6 or 8 mm offset) and fix the coracoid block on the drill guide


As an option the wedged profile plate can be placed between bone and guides.



- Position the glenoid drill guide with the coracoid block on the glenoid rim and drill two guidewires through the guide into the glenoid



- Remove the drill guide and shape the graft with a rongeur, if necessary
- Measure the required screw length using the measuring device
- Use a 2.75 mm drill over the guidewire to predrill the glenoid



- The final fixation is achieved by two cannulated 4.5 mm screws

Ordering Information

Product Description	Item Number
Implants / Disposables	
Wedged profile plate	AR-8111
Low Profile Screw™, Ti, 4.0 x 30 mm, cannulated, long thread, cancellous	AR-8740-30PTL
Low Profile Screw™, Ti, 4.0 x 32 mm, cannulated, long thread, cancellous	AR-8740-32PTL
Low Profile Screw™, Ti, 4.0 x 34 mm, cannulated, long thread, cancellous	AR-8740-34PTL
Low Profile Screw™, Ti, 4.0 x 36 mm, cannulated, long thread, cancellous	AR-8740-36PTL
Low Profile Screw™, Ti, 4.0 x 38 mm, cannulated, long thread, cancellous	AR-8740-38PTL
K-wire trocar tip 1.3 mm, 150 mm length, for 4 mm screws	AR-14513K
K-wire, 1.4 mm x 300 mm, for 4 mm screws	KW02-300-14.A
Screwdriver, T15 hexalobe, cannulated, for 4 mm screws	AR-8943-09
Low Profile Screw™, Ti, 4.5 mm x 30 mm, cannulated, partially threaded	AR-8945-30PT
Low Profile Screw™, Ti, 4.5 mm x 32 mm, cannulated, partially threaded	AR-8945-32PT
Low Profile Screw™, Ti, 4.5 mm x 34 mm, cannulated, partially threaded	AR-8945-34PT
Low Profile Screw™, Ti, 4.5 mm x 36 mm, cannulated, partially threaded	AR-8945-36PT
Low Profile Screw™, Ti, 4.5 mm x 38 mm, cannulated, partially threaded	AR-8945-38PT
Low Profile Screw™, Ti, 4.5 mm x 40 mm, cannulated, partially threaded	AR-8945-40PT
0.062" (1.57 mm) guidewire, 7" (177.8 mm), for 4.5 mm screws	AR-8941-7
0.062" (1.57 mm) guidewire, 12" (304.8 mm), for 4.5 mm screws	AR-8941-12
Screwdriver shaft, cannulated, 3.5 mm hex, for 4.5 mm screws	AR-8100D
Accessory Instruments	
Glenoid bone augmentation case	AR-8100C
Screw caddy for mini open shoulder - glenoid bone augmentation screw caddy	AR-8100C-SC
Swan retractor, right	AR-8102R
Swan retractor, left	AR-8102L
Gelpi retractor	AR-8104
Gelpi retractor for arthroscopic use	AR-8104A
Nerve & fascia retractor	AR-8101
Bended blade retractor, 18 mm	AR-8100-18
Bended blade retractor, 26 mm	AR-8100-26
Parallel drill guide, 6 mm offset	AR-7000-04
Parallel drill guide, 8 mm offset	AR-7000-05
Handle, drill guide, 6.50" (16.5 cm) long	AR-9215-1-01
Screw length gage, glenoid bone loss	AR-7000-06
Coracoid drill guide	AR-7000-07
Fukuda style retractor	AR-7000-08
Drill, 2.75 mm, 0.066" (1.67 mm) cannulation	AR-7000-14
Bio-Tenodesis™ screw drill, 4.0 mm	AR-1204D
Driver handle with AO connection, medium	AR-13421AO
Chondro osteotome (100 mm)	AR-1767
300 sagittal saw blade, angled, 19 x 10 x 0.6 mm	AR-300-450S
Temporary compression device, cannulated	AR-14023TCDC
Optional	
Mini open shoulder retractor, 21 mm	AR-8100-21
Modular soft tissue retractor atraumatic replacement paddle, 50 mm, right	AR-8170-50DR
Modular soft tissue retractor atraumatic replacement paddle, 50 mm, left	AR-8170-50DL
Osteotome blade, glenoid bone loss	AR-7000-01
Osteotome blade shield	AR-7000-02
Dovetail meniscal allograft osteotome handle	AR-2961

Please note that not all products advertised in this brochure/surgical technique guide may be available in all countries. Please ask the Arthrex Customer Service or your local Arthrex Representative before ordering if the desired product is available for delivery. Thank you very much.

Developed in conjunction with Dr. G. Di Giacomo, Dr. A. Costantini, Dr. A. DeVita, Prof. Ph. Hardy, Dr. N. Gravelleau, Dr. J. Barth, Dr. S. Lichtenberg.



This description of technique, including any post-op protocol, is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's Directions For Use.

View U.S. Patent information at www.arthrex.com/corporate/virtual-patent-marking

© Arthrex GmbH, 2016. All rights reserved. | www.arthrex.com | LT2-0550-EN_G