


FlipCutter[®] II

*Arthroscopic Drilling
with the Push of a
Button*

*FlipCutter II is changed
from a straight-tipped drill pin
to a retrograde reamer by
simply pushing the button and
sliding the hub forward*




*The RetroConstruction™ Drill Guide
Instrument Set includes anatomic ACL,
PCL, Meniscus and OATS marking hooks
to facilitate drilling with FlipCutter*



Use FlipCutter II for ...

Anatomic, minimally invasive ACL Reconstruction with RetroConstruction Guides



PCL Reconstruction



Meniscal Root Repair



Ordering Information

FlipCutter II, 6 mm - 13 mm AR-1204AF-60 - 130



For more information go to:

<http://flipcutter.arthrex.com>