

Revolutionizing Small Joint Orthopaedic Surgery with

FiberWire®

FiberLoop® with needle (above)

Mini BioComposite SutureTak®
w/2-0 FiberWire

Small Bone FASTak® Suture Anchor
w/2-0 FiberWire

3.5 mm and 5 mm Corkscrew®
Suture Anchor w/FiberWire

2-0 FiberWire

Micro and Mini Corkscrew® FT
w/2-0 FiberWire

Bio-Corkscrew® FT Suture Anchor
w/#1 FiberWire

3.5 mm Corkscrew® FT Suture
Anchor w/#1 FiberWire

2.5 mm Bio and PEEK
Mini PushLock®

3.7 mm Bio-Corkscrew Suture
Anchor w/#0 FiberWire

FiberWire Scissor, Small



FiberWire


ORTHOPAEDIC COMPOSITE SUTURE

FIBERWIRE IS COMPOSED OF A MULTI-STRANDED LONG CHAIN ULTRA-HIGH MOLECULAR WEIGHT POLYETHYLENE (UHMWPE) CORE WITH A POLYESTER AND UHMWPE BRAIDED JACKET. THIS CONFIGURATION GIVES FIBERWIRE SUTURE EXCEPTIONAL STRENGTH, VIRTUALLY ELIMINATING INTRAOPERATIVE SUTURE BREAKAGE.

Knot Strength Comparison

Peak Load in KG Force (average)





















FiberWire
Ethibond Excel



















See catalog for complete FiberWire offering



FiberWire, FiberLoop and FiberWire Scissor

Suture Size	Suture Length	Description	Catalog Number	Needle Dimensions	Actual Size
2-0 (3 metric)	18"	FiberWire w/Tapered Needle (heavy duty), 1 strand (blue)	AR-7242	26.5 mm 1/2 circle	
2-0 (3 metric)	18"	FiberWire w/Tapered Needle, 1 strand (blue)	AR-7220	17.9 mm 3/8 circle	
2-0 (3 metric)	38"	FiberWire, 1 strand (blue)	AR-7221	-	
3-0 (2 metric)	18"	FiberWire w/Diamond Point Needle, 1 strand (blue)	AR-7225	26.2 mm 3/8 circle	
3-0 (2 metric)	18"	FiberWire w/Tapered Needle, 1 strand (blue)	AR-7227-01	15 mm 3/8 circle	
3-0 (2 metric)	18"	FiberWire w/Reverse Cutting Needle, 1 strand (blue)	AR-7227-02	16.3 mm 3/8 circle	
4-0 (1.5 metric)	18"	FiberWire w/Diamond Point Needle, 1 strand (white)	AR-7248	18.7 mm 3/8 circle	
4-0 (1.5 metric)	18"	FiberWire w/Diamond Point Needle, 1 strand (blue)	AR-7228	18.7 mm 3/8 circle	
4-0 (1.5 metric)	18"	FiberWire w/Tapered Needle, 1 strand (blue)	AR-7230-01	12.3 mm 3/8 circle	
4-0 (1.5 metric)	18"	FiberWire w/Reverse Cutting Needle, 1 strand (blue)	AR-7230-02	11.9 mm 3/8 circle	
0 (3.5 metric)	38"	FiberWire w/Tapered Needle, 1 strand (blue)	AR-7250	22.2 mm 1/2 circle	
0 (3.5 metric)	38"	FiberWire w/Diamond Point Needle, 1 strand (blue)	AR-7251	22.2 mm 1/2 circle	
4-0 (1.5 metric)	6"	FiberLoop w/Tapered Needle, 1 strand (white)	AR-7249-12	12.7 mm 1/2 circle	
4-0 (1.5 metric)	10"	FiberLoop w/Tapered Needle, 1 strand (white)	AR-7249-20	12.7 mm 1/2 circle	
4-0 (1.5 metric)	6"	FiberLoop w/Tapered Needle, 1 strand (blue)	AR-7229-12	17.9 mm 3/8 circle	
4-0 (1.5 metric)	10"	FiberLoop w/Tapered Needle, 1 strand (blue)	AR-7229-20	17.9 mm 3/8 circle	
2-0 (3 metric)	30"	FiberLoop w/Diamond Point Needle, 1 strand (blue)	AR-7232-01	48 mm 1/2 circle	
2-0 (3 metric)	24"	FiberLoop w/Diamond Point Needle, 1 strand (blue)	AR-7232-02	26.2 mm 3/8 circle	
2-0 (3 metric)	13"	FiberLoop w/Diamond Point Straight Needle, 1 strand (blue)	AR-7232-03	64.8 mm	
-	-	FiberWire Scissor, small	AR-11797	-	

Small Joint Suture Anchors with FiberWire & Needles

Suture Size	# of Sutures	Description	Disposable Kit	Cat Number	Anchor Dimensions	Actual Size
2-0 (3 metric)	1	Micro BioComposite SutureTak	AR-1320DSC	AR-1320BCNF	2.4 mm x 6.5 mm	
2-0 (3 metric)	1	Mini BioComposite SutureTak	AR-1322DSC	AR-1322BCNF	2.4 mm x 8.5 mm	
1 (4 metric)	1	BioComposite SutureTak w/Needles	AR-8934DSC	AR-8934BCNF	3 mm x 14 mm	
2-0 (3 metric)	1	Micro Bio-SutureTak	AR-1320DSC	AR-1320BNF	2.4 mm x 6.5 mm	
2-0 (3 metric)	1	Mini Bio-SutureTak	AR-1322DSC	AR-1322BNF	2.4 mm x 8.5 mm	
0 (3.5 metric)	2	Bio-SutureTak w/Needles	AR-8934DSC	AR-8934BNF	3 mm x 14 mm	
2-0 or 0 recommended	0	2.5 mm Bio-PushLock	AR-1322DSC	AR-8825B	2.5 mm x 7 mm	
2-0 or 0 recommended	0	2.5 mm PEEK PushLock	AR-1322DSC	AR-8825P	2.5 mm x 7 mm	
0 (3.5 metric)	2	Bio-Corkscrew w/Needles	Reusable Tap	AR-1923BNF	3.7 mm x 17.9 mm	
1 (4 metric)	2	Bio-Corkscrew FT Suture Anchor	Reusable Tap	AR-8927BNF	5.5 mm x 15 mm	
2-0 (3 metric)	1	Micro Corkscrew FT Suture Anchor	K-wire included	AR-1318FT	2.2 mm x 4 mm	
4-0 (1.5 mm)	1	Micro Corkscrew FT Suture Anchor	K-wire included	AR-1318FT-40	2.2 mm x 4 mm	
2-0 (3 metric)	1	Mini Corkscrew FT Suture Anchor	K-wire included	AR-1319FT	2.7 mm x 7 mm	
0 (3.5 metric)	2	3.5 mm Corkscrew FT Suture Anchor	K-wire included	AR-1915FT	3.5 mm x 10 mm	
2-0 (3 metric)	1	Small Bone FASTak Suture Anchor	(self-drilling)	AR-1322-752SF	2.4 mm x 7.5 mm	
0 (3.5 metric)	2	Corkscrew Suture Anchor	(self-drilling)	AR-1915SNF	3.5 mm x 12 mm	
1 (4 metric)	2	Corkscrew Suture Anchor	(self-drilling)	AR-1920SNF	5 mm x 15 mm	