


# FiberTape<sup>®</sup> Cerclage System

## Fracture Management

The FiberTape Cerclage System is a universal, nonmetallic, suture-based solution for fracture fixation and stabilization. Its design and mechanical properties make FiberTape cerclage an ideal choice for fracture management over metal wires and cables.


- 100% radiolucent – no radiographic interference from cerclage
- No concerns of metal mismatch with adjacent implants
- No sharp ends or broken wires to potentially harm the patient or medical staff
- No metal debris migration or concerns of metallosis

## Trusted biomechanical performance and strong fixation




### References

1. Arthrex, Inc. Data on file (APT 3197). Naples, FL; 2017.
2. Arthrex, Inc. Data on File (APT 04426). Naples, FL; 2019




**Sign up to receive weekly Arthrex What's New emails** for exclusive access to new surgical technique videos, medical education events, animations, and presentations before they appear on our website. Go to <https://arthrex.info/ArthrexLiveTV> or text **WNASM to 278479** to receive videos via mobile alerts.\*

\*Message and data rates may apply when sending and/or receiving text messages. Mobile alerts available only in the United States.



Low-profile tapes lay flat under or above plating

Allows plate to be positioned flush against bone.



Broad footprint compression around bone

Offers the ability for tensionable prophylactic cerclage, which may help decrease the incidence of intraoperative fractures when dealing with poor bone quality.


Tensioner helps achieve the right amount of tension for fracture reduction.

## Ordering Information

Product Description	Item Number
FiberTape® Cerclage Suture	AR-7267
TigerTape™ Cerclage Suture	AR-7267T
Tensioner	AR-7800
Tensioner Ratcheting Handle	AR-7801
Passing Hook, medium	AR-7806
Passing Hook, larger	AR-7807
FiberTape Cerclage Tray	AR-7800C
<b>Cerclage Tensioner Set</b> Includes Tensioner, Tensioner Ratcheting Handle, Medium Passing Hook, Large Passing Hook, and FiberTape Cerclage Tray	AR-7800S
<b>Disposable Passing Needle and FiberLinks</b>	
Passing Needle	AR-7816
FiberLink™ Suture	AR-7235
TigerLink™ Suture	AR-7235T

Streamlined instrumentation fits in a single compact tray

One simple tensioning device and no bulky crimpers or cable cutters.



Visit us at Booths 14 and 15.

[www.arthrex.com](http://www.arthrex.com)

© 2019 Arthrex, Inc. All rights reserved. FL1-000108-en-US\_A

